



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Compresor dentar Energy 24/150

Pret:

4.650,00 lei TVA inclusa



Descriere produs:

- Putere electrica: 1100W / 1,5 CP
- Presiune maxima de lucru: 8 bar
- Presiune normala de lucru: 5-7 +/- 0,5 bar
- Debit de aer: 150l/min
- Capacitate rezervor (butelie): 24l
- Nivel de zgomot: 75 dB
- Dimensiuni: 41 x 40 x 65 cm (latime x adancime x inaltime)
- Masa: 31 Kg
- 1400 rpm
- 1 cilindru
- Potrivit pentru 1 maxim 2 unituri (daca aspiratia este chirurgicala, prin pompa separata)