



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Pungi Sterilizare Pupinel 76X254mm

Pret:

152,00 lei TVA inclusa



Descriere produs:

Pungi sterilizare pupinel, autosigilante, dimensiuni 76x254mm

Prezentare:

- Pungi sterilizare pupinel 76x254mm - 1 x 100 pungi

Mod de folosire: vezi sectiunea download