



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Ketac Universal Aplicap



Pret:

679,99 lei TVA inclusa

Descriere produs:

Ketac Universal Aplicap - ciment ionomer de sticla pentru restaurari rapide, usoare si rezistente

Avantaje:

- Aplicare intr-un singur pas - nu necesita conditioner sau strat protector (lac, varnish)
- Rezistenta la uzura - ideal pentru restaurarile de clasa I, II, V
- Elibereza fluor continuu, timp de 24 luni
- Se aplica bulk - nu necesita stratificare
- Adaptare marginala usoara
- Autoadeziv
- Autopolimerizabil
- Disponibil in 6 nuante: White (alb), A1, A2, A3, A3.5, A4
- Radioopac
- Manevrabilitate usoara si aplicare rapida
- Adera la cavitate in timpul aplicarii, dar nu adera la instrumentele dentare

Indicatii:

- Ca liner inainte de aplicarea obturatiilor de compozit



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

- Refaceri de bont
- Restaurari pe dintii temporari
- Restaurari de Clasa I si Clasa a II-a*
- Restaurari de Clasa a V-a
- Restaurari temporare
- Obturatii minim-invazive

Prezentare:

- Ketac Universal Aplicap - 1 cutie x 50 capsule