



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Affinis Putty Soft



Pret:

565,00 lei TVA inclusa

Descriere produs:

Silicon de aditie tixotrop, excelenta afinitate in mediul umed datorita adaosului de surfactanti, fluiditate crescuta

Indicatii

- acuratete in reproducerea campului protetic, astfel incat amprenta nu va prezenta lipsuri sau distorsiuni ale marginilor sale

Prezentare

- 300 ml baza
- 300 ml catalizator