



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Snur Retracto impregnat

Pret:

68,34 lei TVA inclusa



Descriere produs:

Fir pentru retractie gingivala impregnat cu clorura de aluminiu.

Indicatie

- retractie gingivala
- mici hemoragii gingivale

Prezentare

- 1 cutie x 1,83m
- 4 dimensiuni: X-fine, Fine, Mediu si Gros
- impregnat