



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehncaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

---

## Lampa examinare chirurgicala Micare JD1700 podea



**Pret:**

**5.950,00 lei TVA inclusa**

### **Descriere produs:**

Lampa chirurgicala mobila, cu LED pentru chirurgie orala si maxilo-faciala, chirurgie plastica/estetica.

Montare pe podea.

Design modern si practic.

### **Caracteristici:**

Iluminare cu LED, lumina rece, luminozitate ridicata

Iluminare omnidirectionala

Comutator tactic inteligent

Consum redus de energie

Intensitatea luminii reglabila (5200-80000 lux)



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

---

Temperatura culorii: 4000-6500K

Tensiune intrare: 100-200V ~ 50 / 60Hz 12V / 2A-2.5A

Tensiune ieșire: 12V / 2A ~ 2.5A

Putere: 24W

Durata de viata a becului: 50000 ore

Dimensiunea suportului: 37 cm, inaltime 8.5 cm