



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

---

## Micromotor MB1-A

**Pret:**

**489,99 lei TVA inclusa**



### **Descriere produs:**

Micromotor pneumatic, conexiune Midwest (4 orificii).

Prevazut cu tub si conector pentru irigare externa.

Viteza: 22.000 rpm la o presiune de alimentare cu aer de 2.5 bar (245 kPa).

Consum maxim de aer: 42 l/min.

Reglaj turatie si inainte/inapoi.

Necesitate de lubrefiere: zilnica

Autoclavabil la 134C.

Procedura de autoclavare:



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

---

- stergeti impuritatile cu servetel umed

- lubrefiati

-introduceti in punga de sterilizare

- autoclavati la 121C sau 134C