



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Aparat detartraj UDS-P LED

Pret:

1.065,00 lei TVA inclusa



Descriere produs:

Aparat de detartraj cu ultrasunete

Caracteristici

Functionalitate dependenta fata de unit, necesita o conectare la unitul dentar

Control al apei

Piesa de mana din titan, detasabila si autoclavabila

Furtun confectionat din gel siliconic, moale si durabil.

Control digital

Intensitatea frecventei se regleaza automat in functie de situatia data

Putere: 3-20W

Frecventa: 28KHz



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Funcții: detartraj, perio, endo

Pachetul contine 6 anse si 4 ace endo din titan