



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Gingifast rigid 2x50ml

Pret:

218,90 lei TVA inclusa



Descriere produs:

Silicon de aditie de laborator indicat pentru masca gingivala.

Avantaje:

- Usor de ingrijit permite o definire de acuratete a marginilor
- Indicat pentru protetica pe implant
- Precizie si stabilitate dimensionala

Caracteristici:

- Silicon de aditie ambalat in cartuse

Prezentare:

Cartus 2 x 50ml + 10ml separator + 12 canule de amestec galbene + 12 varfuri intraorale, galbene.