



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

---

## Gingifast elastic 2x50ml

**Pret:**

**220,60 lei TVA inclusa**



### **Descriere produs:**

Silicon de aditie de laborator pentru masca gingivala (elastic).

### **Avantaje:**

- Modulul elastic permite indepartarea si reinsertia
- Rezultate estetice excelente datorita proprietatilor de transluciditate
- Precizie si stabilitate dimensionala
- Folosit in tehnica directa si indirecta

### **Caracteristici:**

- Silicon de aditie ambalat in cartuse

### **Prezentare:**

- Gingifast Elastic (2 x 50 ml), 10ml separator, 12 canule de amestec galbene, 12 varfuri intraorale galbene