



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

CUPTOR CALCINARE STC

Pret:

6.000,00 lei TVA inclusa



Descriere produs:

Temperatura maxima: 1100° C

Carcasa: inox

Capacitate 4 litri

Controler digital

Precizie stabilizare: $\pm 2^\circ$ C



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Alarma cu timer

Izolatie termica cu fibra ceramica

Rezistente din KANTHAL-Suedia

Protectie la termocuplu ars

Dimensiuni interne 155x210x110mm

Putere consumata: 1.4 KW

Alimentare:230 V