

## HERAENIUM P

**Pret:**

**1,50 lei TVA inclusa**



### **Descriere produs:**

**Compoziție:** Co 59%, Cr 25%, W 10%, Mo 4 %, Si 1%, Mn 0,8%

### **Caracteristici:**

- Densitate 8,8g/cm<sup>3</sup>
- Coeficient de dilatare termica 13,8μm/mK (25 – 500 grade C)
- Duritate 330 HV
- Interval de topire: 1305 grade C-1400 grade C
- Temperatura de turnare 1550 grade C.

### **Prezentare:**

Gram.