



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

## CERAMIC MASTER PRESS VOP

Pret:

**34.500,00 lei TVA inclusa**



### Descriere produs:

#### Caracteristici:

- + Destinat arderii tuturor tipurilor aliaje dentare si ceramica pura si presata
- + nu este necesar compresorul pentru procesul de presare
- + 100 de programe individuale de operare
- + Proces controlat in totalitate electronic
- + Sistem imbunatatit pentru autodiagnosticare si autocalibrare a temperaturii
- + Pompa de vacuum inclusa PVM122

Consum: 1700W (cu pompa de vacuum) / 1550W (fara pompa de vacuum)

Putere de absorbtie a pompei - 22l/min

Vacuum ajustabil 0.1 - 0.92 bar



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

---

Temperatura maxima - 1200 °C

Dimensiuni camera de ardere -  $\Phi$  95 mm, H 60 mm

Dimensiuni (L/l/h) 370x380x740

Masa: 30 kg