



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Ciment dual Bifix QM 10g

Pret:

645,00 lei TVA inclusa



Descriere produs:

Sistem compozit de cimentare universal cu polimerizare duala.

Indicatii:

- Fixarea aluziunilor pe baza de ceramici, in crustele de tip Inlay Onlay ceramicele puntilor puntilor Maryland*

Avantaje:

- Polimerizare duala Bifix QM (cu Ceramicele Bifix) permite o aplicare simpla si rapida la portelan*

Prezentare:

- Bifix QM Universal syringa 10g*