



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

---

## Variotime Dynamix Heavy Tray 2 x 380 ml

**Pret:**

**420,00 lei TVA inclusa**



### **Descriere produs:**

Silicon de tip A, hidroactiv, pe baza de vinilpolisiloxan ce ofera precizie pentru toate tipurile de tehnici de amprentare, consistenta medie la mixare.

### **Indicatii:**

- Tehnica de amprentare cu o singura vascozitate - materialul de amprenta unicomponent se va folosi simultan atat in lingura cat si in jurul bontului.

### **Avantaje:**

- Timp de lucru variabil: intre 1' si 2'30";
- Timp de priza in cavitatea bucala constant: 2'30".

### **Caracteristici:**

- Deformare scazuta a amprentei;
- Tixotropie;



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

---

- Reproducerea detaliilor de suprafata cu acuratete.

### **Mod de prezentare**

1 x cutie Variotime Heavy Tray 2 cartuse x 380 ml