

Alginat Ruthinium



Pret:

31,00 lei TVA inclusa

Descriere produs:

Caracteristici

- Alginat cromatic de inalta precizie cu priza normala, marcat cromatic.
- Stabilitate dimensionala excelenta a amprentei timp de 5 zile, daca este depozitata corect, in container inchis ermetic.
- Consistenta supla.
- Material non toxic.
- Timp de priza: 2 min si 20 sec cu apa la temperatura de 23 C.
- Preparare usoara

Indicatii

- alginat pentru amprentare in scop ortodontic

Prezentare



Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

- punga: 453g