



## Aliaj dentar Niadur



**Pret:**

**1,05 lei TVA inclusa**

### **Descriere produs:**

Aliaj Ni-Cr tip 3, pentru scheletul metalic al lucrarilor de ceramica. Grosimea peretelui la o singura coroana: min. 0.3mm, la o punte: min. 0.4mm.

### **Caracteristici tehnice:**

Compozitie Ni:62%, Cr 24,50%, Mo 10,30%, Si 1,6-2%, Fe 1,5%.

Coeficient de expansiune(25-500 C) 14.0  $\mu\text{m}$ .

Densitate 8.2g/cm<sup>3</sup>.

Interval topire: 13000C.

Duritate Vickers 200 HV10.

Temperatura de turnare 15000C.

Temperatura de preincalzire: 900-9500C.

Alungire la rupere: 4.8%.



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

---

Limita de intindere: 380 N/mm<sup>2</sup>.

Modul elasticitate: 200.000N/mm<sup>2</sup>.

Usor de prelucrat si lustruit.

Pret / gram