



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

---

## Seringi cu ac 2 ml (100 buc)

**Pret:**

**43,49 lei TVA inclusa**



### Descriere produs:

- seringi cu ac pentru clatire.
- 100 buc/cutie
- **2 ml cu ac 23G x 1"**

Seringile si acele de seringa marca Demotek sunt fabricate din materiale (polipropilena si cauciu natural) pentru uz medical de inalta calitate, in concordanta cu normele BS EN ISO 9002 si detin marca CE conform Directivei Europene 93/42/EEC.

Producator: Demophorius Limited Cipru, UE

Seringi din 3 componente (piston prevazut cu garnitura din cauciu natural pentru uz medical)

Dimensiuni seringi: 2ml

Dimensiuni ace: 2ml - 23G

Disponibile in varianta luer-slip

Metoda de sterilizare: ETO (Ethylene Oxide)



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

---