



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Kit HAWE posterior restauration

Pret:

690,56 lei TVA inclusa



Descriere produs:

Hawe posterior restoration kit contine unele dintre cele mai utilizate produse HAWE - penele SICAMORE, Optidisc, Optishine, Optistrip, si unele dintre cele mai inovative produse de pe piata - Metafix si Opti1Step Polisher.

Continut

- 20 x Metafix Small
- 30 x Metafix Medium
- 20 x Metafix Large
- 50 x Pene interdentare Turquoise
- 50 x Pene interdentare Orange
- 50 x Pene interdentare Albe
- 3 x Periuta lustruit Optishine
- 2 x Guma lustruit Opti 1 Step tip flacara
- 2 x Guma lustruit Opti 1 Step tip cupa
- 50 x Banda finisare OptiStrip 4 mm
- 50 x Banda lustruire OptiStrip 4 mm
- 2 x Mandrina OptiDisc Standard



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

- 2 x Mandrina OptiDisc scurta
- 20 x Optidisc extradur 12,6 mm
- 20 x Optidisc duritate medie 12,6 mm
- 20 x Optidisc fine 12,6 mm
- 20 x Optidisc extra-fine 12,6 mm