



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehncaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Aparat de detartraj incorporabil UDS-N3 LED

Pret:

920,00 lei TVA inclusa



Descriere produs:

Caracteristici

- incorporabil in unit
- urmarire automata a frecventei
- piesa detasabila cu lumina LED

UDS-N3 LED - aparat detartraj cu ultrasunete pentru unit

Avantaje:

- Piesa de mana detasabila cu LED, autoclavabila
- Urmarire automata de frecventa
- Functii: detartraj, endo, perio
- Varfuri compatibile incluse: 1 x E1T, 2 x G1T, 1 x G2T, 1 x G4T, 1 x P1T

Date tehnice:

- Tensiune unitate principala: 24V, 50-60Hz, 1.3A
- Consum: 3-20W



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

- Frecventa: 28kHz \pm 3kHz
- Presiune debit apa: 0.1-5 bari (0.01-05MPa)
- Greutate unitate principala: 0.2kg
- Dimensiuni: 74 x 56 x 38mm

Pachetul contine

- UDS-N3 LED - 1 aparat
- 6 anse incluse: 1 x E1T, 2 x G1T, 1 x G2T, 1 x G4T, 1 x P1T