



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro  
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

## Endo motor Pro Wireless

Pret:

**2.350,00 lei TVA inclusa**



### Descriere produs:

Endo Motor cu Apex locator EDO-1 Pro Wireless.

Endo Motor cu Apex locator EDO-1 Pro Wireless

EDO-1 Pro Wireless

Tip Portabil Memorie 10 Setari memorate Rotatii In sensul acelor de cesornic/Anti-Clockwise/Reciproc  
Torque Control Auto Reverse/Auto Stop Baterie Reincarcabila/Inlocuibila Incarcare Charging Base/Micro  
USB Display OLED Blue Tooth Available Contra Angle 1:1 Apex Locator Control Motor Working Mode  
Apex locator de precizie superioara necesar in detectarea apexului radicular.

Apex locator, necesar de precizie superioara necesar in detectarea apexului radicular.

Apex locator, necesar de precizie superioara necesar in detectarea apexului radicular.

Player video