



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Surgitip 2.5mm

Pret:

245,00 lei TVA inclusa



Descriere produs:

Aspiratie si precizie in timpul procedurilor microchirurgicale

Canula ofera o priza buna si este usor de manevrat, facand posibila aspiratia tintita si in mediile operatorii dentare restrictionate

Steril, ambalat individual, gata pentru utilizare imediata.

Canule de aspiratie chirurgicala - 2.5 mm.

Confortabile pentru pacient.

Au o inalta calitate si pot fi sterilizabile in autoclav cls. B.

20 buc/cutie + 1 adaptor