



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Mazic D Temp

Pret:

1.350,00 lei TVA inclusa



Descriere produs:

Rasina tip acrilat fotopolimeric, pentru imprimarea coroanelor si puntilor.

Compatibila cu imprimantele 3D de intarire cu lumina UV.

Caracteristici si beneficii

Estetica dentara naturala

Echilibru intre opacitate si translucenta coroanelor si a puntilor.

Rezistenta ridicata

Precizie mare

Viteza de imprimare

Solubilitate reduca in apa si absorbtie.

Culoare: verde opac inainte de imprimare



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Vascozitate: 1100~1300 cP

Mod de prezentare

500 g/1 Kg Mazic D Temp