



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Varfuri intraorale compatibile tip 1

Pret:

44,99 lei TVA inclusa



Descriere produs:

- Pentru Rebuilda DC in seringa Quick Mix, Bifix SE, Bifix QM. Se potrivesc cu varfurile de mixare HP06.
- 50 buc/punga
-