



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Tavita rasina pentru imprimanta Accufab L4D Shining 3D

Pret:

650,00 lei TVA inclusa



Descriere produs:

Tank / Tavita rasina pentru imprimanta Accufab L4D de la Shining 3D

Avantaje:

Durabil si solid

Rama metalica, ofera rezistenta si o durata de viata extinsa

Capac reutilizabil ce permite inmagazinarea usoara a rasinii

Pentru AccuFab-L4D

Prezentare

Tavita rasina L4D - 1 bucata



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Capac Tavita – 1 bucata