



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Max Filltemp 28g



Pret:

50,16 lei TVA inclusa

Descriere produs:

Ciment de umplere temporara pe baza de rasina sintetica.

Ideal pentru etansarea inlay-urilor si coroanelor, adaptandu-se perfect.

Nu necesita amestecare.

Usor de manipulat datorita consistentei sale mai fluide si mai eficient datorita duritatii dupa prizare.

Prezentare

1 x flacon 28g