



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Pungi autosigilante 83x159mm

Pret:

40,00 lei TVA inclusa



Descriere produs:

Pungi autosigilante din polietilena si hartie. Panoul indicator de culoare se schimba cand a avut loc sterilizarea.

Dimensiune: 83mm x 159mm

Prezentare

200 bucati