



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Dezinfectant concentrat Eurosept Xtra



Pret:

194,00 lei TVA inclusa

Descriere produs:

Concentrat pentru dezinfectia si curatarea manuala a instrumentelor generale si chirurgicale intr-o baie de imersie. (pH >10) Potrivit pentru utilizare intr-o baie cu ultrasunete.

Timpe de aplicare si dilutie:

Bacteriile conform EN 13727, EN 14561: 1,5% (15 ml pe litru) 15 minute, 1,0% (10 ml pe litru) 60 minute;

Candida albicans conform EN 13624, EN 14562: 1,5% 15 minute, 1,0% 60 minute;

Bacteriile si Candida albicans intr-o baie cu ultrasunete conform EN 14561 si EN 14562: 1,0% 5 minute;

Virucid limitat (virusuri incapsulate, cum ar fi HBV, HIV, HCV) EN 14476 si EN 17111: 1,5% 15 minute;

TBC si micobacterii (M. terrae si M. avium) conform EN 14348 si EN 14563: 4,0% (40ml pe litru) 30 minute, 3,0% (30ml pe litru) 60 minute.

Prezentare



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Flacon 2L