



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Ansa endo Sonic File Holder E2

Pret:

90,00 lei TVA inclusa



Descriere produs:

Sonic File Holder E2 95° este utilizat in tratamentul dintilor anteriori si canalelor meziale posterioare. Varf compatibil cu dispozitivele cu ultrasunete EMS / WOODPECKER / MECTRON.

Autoclavare cu abur la 134 de grade. Durata recomandata a sterilizarii: 35-40 minute.

Prezentare

Ansa endo Sonic File Holder E2 95°