



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Thixoflex M

Pret:

52,50 lei TVA inclusa



Descriere produs:

- Silicon de condensare cu vascozitate scazuta. Consistenta medie. Fluid hidrocompatibil.
- Culoare: albastru deschis.
- Tixotrop (caracteristica ce imbunatateste fluiditatea la presiune usoara)

Indicatii

- ideal in tehnica simplei amprente cu doua vascozitati
- tehnica dublei amprente
- mucostatic (necompresiv pe tesaturile moi)
- fluiditate ridicata in timpul lucrului
- timpi de pregatire si lucru variabili in functie de cantitatea de Indurent adaugata

Prezentare

- 140 ml (tub)



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273
