



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehncaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Kit Technique 12 cleme diga K-CLAMP Shinhung Coreea

Pret:

990,00 lei TVA inclusa



K-CLAMP TECHNIQUE SET

SET Includes: 12 clamps, 1 punch,
1 forcep, 1 frame

Descriere produs:

Sistem de diga K-clamp: 12 cleme mate, rama diga metalica, cleste aplicare diga, cleste perforator diga.

Materialul metalic este mat non-reflectorizant, ceea ce ofera o vizibilitate sporita in timpul procedurilor dentare, reducand reflexiile luminoase.

Design ergonomic si eficient, pentru a va oferi cel mai bun acces si vizibilitate a locului izolat. Elasticitatea face K-Clamps extrem de durabile si usor de manevrat. Clemele ofera un confort superior pentru pacient si previne traumatismele nedorite ale tesutului moale din jur.

Otel inoxidabil fabricat pentru a rezista la coroziune in conditii intraorale.

AUTOCLAVABILE la 121°C

Coduri cleme:

201, 202 , 203, 204, 205, 206, 207, 208, 209, 12, 210, 211



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Prezentare

Set de diga x 15 piese