



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Amaris Seringa 4g TL EXP 28.09.2024

Pret:

189,00 lei TVA inclusa



Descriere produs:

Indicatii

- Restaurari estetice ale dintilor frontali clasele III, IV si V
- Restaurari estetice ale dintilor posteriori clasele I, II si V
- Restaurari ale dintilor frontali cu leziuni de origine traumatica
- Fatetari directe si indirecte in compozit, fatatari pe dintii frontali colorati
- Corectari de forme si culori cu scop estetic
- Reparatii de proteze in ceramica estetica si fatetari

Avantaje

- Lustruire simpla, rezistenta la abraziune si o stralucire indelungata
- Consistenta ne-adeziva, rezistenta durabila la lumina de operatie
- Doua etape simple, putin material necesar
- O stratificare naturala
- Dinamica luminoasa excelenta, fluorescenta naturala
- Stabilitate crescuta a culorii



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

- Acopera tot spectrul culorilor prin combinatia celor 8 componente
-

Mod de prezentare:

- 1 x seringa de 4 g in culorile O1, O2, O3, O4, O5, TL (translucent light), TN (translucent normal), TD (translucent dark)
-