



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Rasina Ghid Chirurgical 1L EXP 10.08.2024



Pret:

367,50 lei TVA inclusa

Descriere produs:

Avantaje

Pentru printarea ghidurilor dentare chirurgicale

Biocompatibila

SG01 - polimer rigid cu flexiune minima, care ii permite sa isi mentina forma in timpul chirurgiilor dentare

Transparent

Proprietati fizice (stare lichida)

Culoare: transparent

Vascozitate (25°C): 527cps

Lungime de unda: 405nm



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Proprietati mecanice (dupa polimerizare)

Duritate Shore D: 84

Modul elastic: 1873MPa

Rezistenta la tractiune: 43MPa

Modul de flexiune: 2047MPa

Rezistenta la inconvoiere: 81MPa

Post-procesare

Timp de curatare: 6 minute

Timp de polimerizare recomandat: 5 minute pentru fiecare fateta

Prezentare

Rasina ghid chirurgical - Flacon 1L