



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehncaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Elite HD Light Body tub

Pret:

209,00 lei TVA inclusa



Descriere produs:

Descriere

Silicon de aditie cu vascozitate mare.

Indicatii

Silicon de aditie cu priza rapida pentru inregistrarea amprentei protezelor fixe.
Tehnica indicata - cea in doua etape, cu linguri de amprenta standard.

Prezentare

- 2 tuburi x90ml.