



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Rasina printare Splint Soft 1L

Pret:

1.725,00 lei TVA inclusa



Descriere produs:

Avantaje

Rezistenta termica si duritate ridicate - permit modelelor polimerizate sa reziste la abraziune si mentin precizia in utilizari multiple

Ideala pentru printarea pentru fabricarea alignere-lor si a altor aplicatii ortodontice

Proprietati fizice (stare lichida)

Culoare: Transparent

Vascozitate (25°C): 1050cps

Lungime de unda: 405nm

Proprietati mecanice (dupa polimerizare)



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Duritate Shore D: 82

Modul elastic: 2003MPa

Rezistenta la tractiune: 1,5 %: >23MPa

Modul de flexiune: >200MPa

Rezistenta la inconvoiere: >5MPa

Post-procesare

Timp de curatare cu etanol 95%. Curatare manuala: 20 secunde

Timp de polimerizare recomandat Fabcure 1/1S: 10 minute | Fabcure 2: 5 minute

Prezentare

Rasina splint ortodontic - Flacon 1L