



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Ciment fotopolimerizabil NX3 WHITE seringa 1.8g EXP 21.11.2024

Pret:

144,00 lei TVA inclusa



Descriere produs:

Ciment rasinic de generatia a treia, cu adeziune excelenta la dentina, smalt, blocuri CAD/CAM, ceramica, portelan, rasina si metal

- fotopolimerizabil, indicat si in cimentarea fatetelor si pivotilor
- compatibilitate atat cu sistemul self-etch cat si cel total-etch

Prezentare

- seringa 1.8g