



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Primsole Putty Soft 2x250ml

Pret:

190,01 lei TVA inclusa



Descriere produs:

Materiale de amprenta (VPS) 2x250ml

AVANTAJE

Rezistent: rezistenta mare la rupere

Reproducere precisa a detaliilor

Stabilitate dimensionala indelungata

Recuperare elastica excelenta

Sigur

Termen de valabilitate 3 ani



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

CARACTERISTICI

Consistentă moale

Chit tip 0 – ISO 4823

Timp îndelungat în cavitate

Duritate finală 65 Shore A

Biocompatibil

Culoare: portocaliu

Prezentare

Primsole Putty Soft - 2 x 250ml