



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Primsole Light Body Fast 2x50ml

Pret:

140,00 lei TVA inclusa



Descriere produs:

Materiale de amprenta (VPS) 2 x 50ml

AVANTAJE

Rezistent: rezistenta ridicata la rupere

Reproducere inalta a detaliilor

Stabilitate dimensionala lunga

Recuperare elastica excelenta

Siguranta

Termen de valabilitate de 3 ani

Combinat cu varfuri de amestecare avansate - usor de blocat



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

CARACTERISTICI

Consistentă tixotropica

Corp usor tip 3 - ISO 4823

Hidrocompatibil

Timp rapid in cavitate

Biocompatibil

Scanabil

Culoare: dark blue

Prezentare

Primsole Light Body Fast - 2 x 50ml, 12 varfuri