



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Lampa foto Eighteeth Curing Pen



Pret:

1.198,00 lei TVA inclusa

Descriere produs:

Trei optiuni de fotopolimerizare: Ramp, Pulse, Mod detectare carii

Mod detectare carii: 600 mW/cm²

Putere de lucru pentru fotopolimerizare: 1000 mw/cm², 1500 mw/cm², 2300 mw/cm²

4 LEDuri polimerizare - 3 Albastre, 1 Violet

Lungime de unda: 385-515 nm

Cap cu rotire completa la 360 grade

Display digital mare

Putere constanta de emisie

Foarte usoara: doar 109 grame!



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Design ergonomic: 24*21*230(mm)

Baterie Li-Ion puternica, rezistenta la >1200 polimerizari

Instructiuni utilizare