



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Clip Flow seringa 1.8gr

Pret:

97,01 lei TVA inclusa



Descriere produs:

Compozit fluid fotopolimerizabil pentru restaurari provizorii

Indicatii:

- Restaurari provizorii ale dintilor anteriori si posteriori si ale coroanelor/antelajului radicular*

Avantaje:

- Rezistenta la presiuni si la palpatare*

Prezentare:

- 1 x seringa x 1.8g*