



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Sigilator Seal 120

Pret:

1.350,00 lei TVA inclusa



Descriere produs:

Aparat pentru sigilarea pungilor termosigilante de ultima generatie.

Design compact si inovator

Necesita numai 3 minute timp de preincalzire

Putere: 100W

Latimea benzii de sigilare (amprenta): 12 mm

Lungime sigilare: 300 mm (dimensiunea maxima a roletelor de pungi termosigilante)

Greutate: 4.8 kg.

Temperatura de sigilare $200\pm 100^{\circ}\text{C}$

Dimensiuni: 500 x400x210 mm



Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@tehnicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273
