

ISODENT LICHID (izolant) 250ml EXP 30 01 2025



Pret:

5,00 lei TVA inclusa

Descriere produs:

Isodent - solutie alginata separare gips

Avantaje:

Aplicare usoara

Uscare rapida

Lichid vascos

Culoare rosu spre rosu deschis

Permite izolarea perfecta a modelelor cu doua straturi de gips precum si a matritelor cu 2 straturi (gips/gips sau gips/rasina)

Nota:



Intrarea Sulfinei 94 Magurele, Ilfov. Tel: 021 312 56 60, E-mail: office@tehnicaldent.ro IBAN: RO17 RNCB
0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Se recomanda inmagazinarea intr-un loc uscat si intunecos, la o temperatura ambientala de 5 – 25°C, in ambalaj bine inchis.

Prezentare:

Solutie alginata separare gips Isodent - 1 flacon x 250ml