




Domnita Ancuta nr 8, Sector 1, Bucharest Tel./Fax: 0040 01 312 56 60, E-mail: office@technicaldent.ro
IBAN: RO17 RNCB 0068 0043 7192 0001 Nr. Inreg. Reg. Com.: J/09/137/99, Cod fiscal: RO 117 58 273

Compresor dentar Dental Oilless 50/230



Pret:

12.389,99 lei TVA inclusa

Descriere produs:

- Tip dental complet carcasat
- Putere 1,5 kW / 2 CP
- Volum butelie: 50l
- Debit aer la iesire: 230 l/min
- Nivel de zgomot 62 dB
- Presiune maxima de lucru 8 bar
- Dimensiuni: 50 x 61 x 102 cm (latime x adancime x inaltime)
- Masa: 82 Kg
- 1400 rpm
- 2 cilindri in V
- Potrivit pentru 2 maxim 3 unituri (daca aspiratia este chirurgicala, prin pompa separata)